


INDIA
Personal Protective Equipment (PPE)
Products



1. Disposable Mask
2. 3 Ply Reusable Mask
3. Cover All (Non Surgical/Surgical)

Disposable Mask (N95 Equivalent)



Content : 100% Polpropylene

Material : Non Woven - Laminated

Color : White

Size : L 6.5 x W 4 Inchs

Moq : 50K Pcs

Finishing: 99.9% Bacterial filtration Efficiency

Price : Refer Price Matrix

Nonsurgical Mask-Nonwoven -Disposable



2 Ply Reusable Mask (100 wash)



Content : 100% Polyester

Material : Woven - Plain Weave

Gsm : 70

Color : White

Size : W 7 3/4 x H 4 Inchs


Moq : 10K Pcs

Price : Refer Price Matrix

Note: Product Compliant for FDA approval

Specification Reusable Mask

Specifications




3-Layer Filtration



NOT INTENDED OR CERTIFIED AS A MEDICAL DEVICE

Key Features Reusable Mask



Style MD Reusable Mask



LAYER 3 : KNIT FABRIC FOR
COMFORT AND SAFETY

COMFORTABLE ELASTIC SUPPORT

FUSING STICKER

LAYER 2 : MELT BLOWN FOR BACTERIAL PROTECTION

LAYER 1 : KNIT FABRIC FOR
COMFORT AND SAFETY

Style ETN Reusable Mask

Powered by  [®]
world technology



Breathable



Washable



Soft Inner Knit fabric



Reusable



Anti Dust Pollution



Sterilized for use



Antimicrobial
N9 Pure silver
Technology



Easy Carry on
ziplock pouch




Melt Blown Layer
for enhanced Bacterial
Protection

Size Chart Reusable Mask

FOR ADULT MEN & WOMEN

SIZE MEASUREMENTS (CM)	MEDIUM	LARGE
Width	21.5	23.5
Side Height	5	6
Height	9.5	10



Multi Pack Reusable Mask



3-Layer Filtration



ALLERGEN



DUST



POLLEN



FUR



CHEMICAL HAZARDS



SPRAY

Fabric Assortment Reusable Mask

DOT ALL OVER PRINT



PLAIN FABRIC



FLORAL ALL OVER PRINT



STRIPE'S PRINT



PLAIN FABRIC



FLOWER ALL OVER PRINT



Fabric Assortment Reusable Mask



Features Reusable Mask



WITH MELT BLOWN
PROTECTION



95% BACTERIA
FILTRATION EFFICIENCY
FOR >3 MICRON SIZED BACTERIA




DUST PROTECTION



BREATHABLE



SOFT ELASTIC LOOPS
FOR EXTRA COMFORT



REUSABLE (UPTO
30 GENTLE WASHES)

Product USP

USP OF BASIC CLOTH MASK

- Fashionable masks
- Anti Pollution
- Anti Dust
- Soft Inner knit fabric for added comfort
- Soft Elastic loops for added comfort
- Reusable
- Washable
- Breathable
- Sterilize for your use

Fabric - knitted soft cotton fabric

ADD ON USP OF AB MASK

- Melt Blown filter layer for enhanced bacterial protection

Maintenance Reusable Mask

- Wash your hands before touching your mask.
- When you put on the mask make sure it covers your nose and mouth and fits comfortably around your face.
 - While you have the mask on in public places, avoid touching it and replace it if it feels damp.
 - Do Not wear the mask for more than 6 hours continuously.
- To remove the mask, remove it from behind, do not touch the mask; store away separately in a closed bin.
 - Do not touch the front side of the mask.
 - Wash mask with soap and water / antiseptic liquid after every use.
 - Wash hands carefully with soap after removing mask.



Washable



Easy
maintenance



Breathable



Soft elastic loops
for extra comfort



WARNING

For individual use only, do not share.

Mask are effective only when used in combination with frequent hand
cleaning with alcohol based hand rub or soap & water.

Non Surgical Gown - Level 1/Level 2/Level 3




Note: Product Compliant for FDA approval

Surgical Gown - Non Woven - Disposable



Note: Product Compliant for FDA approval



Man, Machine & Fashion Blended Melange